
# Need More Than One Software Program?



#### Suites are Created by Selecting Different Product Families

• For each product family, choose the software version you need (Basic, Reg, Pro, Arabidopsis or Density)

#### Advantageous Discounts Offered

• The larger software Suite you buy, the bigger the discount

#### High Quality Hardware Recommended Depending on the Software Suite

- Calibrated optical scanners with special lighting system for accurate scientific measurements of thick samples such as roots, seeds, and needles
- Several software programs of a Suite can share the same scanner
- Positioning systems with trays, and wood core holders for easy and efficient scanning
- Blue reflectance target for wood density analysis
- Portable imaging and calibration Kit available for leaf image acquisition with a camera (camera not included)
- · Self-leveling mount with compact DSLR camera and calibrated fisheye lens for canopy analysis





# **REGENT** Image Analysis for Plant Sciences

## New 2022 Software Programs Available Now!

#### WinRHIZO™

Analyser of Washed Roots and Arabidopsis Seedlings





#### WinFOLIA™

Leaf Area Meter and Morphological Analyser





### WinSCANOPY™

Analyser of Canopy Structure and Solar Radiation








#### WinDENDRO™

Tree-Ring and Wood Density Analyser





#### WinRHIZO™ Tron

Interactive Analyser of Images of Roots Growing in Soil






#### WinSEEDLE™

Seed and Needle Morphological Analyser and Counter






### WinCAM NDVI™

Color Area Meter, NDVI & Basic Morphological Analyser





#### WinCELL™

Analyser of Wood-Cell Structure Parameters





## At Least 5 Good Reasons for Choosing our Products:

## Perpetual Software Licences

- No need to budget licence renewals each year
- Use our products as long as you wish
- No activation and no internet connection required
- 32-bit and 64-bit versions included with each licence

## Frequent updates

- Stay up to date with new trends and features
- Update or upgrade software programs anytime
- Available even for old licences bought in 2000

## Free Technical Support

- Prompt since 1991
- Available via email
- Provided by our competent engineers

## High Quality Hardware

- Optical scanners calibrated for accurate scientific measurements
- Special lighting system for thick samples, e.g. roots, seeds, needles
- Positioning systems with trays and wood core holders for easy and rapid scanning of samples
- Blue reflectance target for wood density analysis
- Several software can share the same calibrated scanner and accessories
- Self-leveling mount with compact DSLR camera and calibrated fisheve lens for canopy analysis



#### **New Advantageous Discounts**

 The more software programs you buy, the bigger the discount